CONTRAfuran™

SENSOfluran CONNECT System

Instructions for Use

Instruction for Use - English	1-12
Mode d'emploi - Français	13-24
Gebrauchsanweisung - Deutsch	25-36
Istruzioni per l'uso - Italiano	37-48
Instrucciones de uso - Español	49-60
Instruções de utilização - Portuguese	61-72
Gebruiksaanwijzing - Dutch	73-84
Käyttöohje - Finnish	85-96
Bruksanvisning - Norsk	97-108
Bruksanvisning - Swedish	109-120
Mode d'emploi - French Canadian	121-132

Instruction for Use

These instructions for use contain information about the **CONTRAfuran** / SENSOfluran CONNECT System, consisting of the **CONTRAfuran** Anesthetic gas canister, the SENSOfluran CONNECT Fill level control unit and associated components.

Please read the instructions for use carefully before performing installation and operation of the **CONTRAfuran / SENSOfluran CONNECT Anesthetic Gas Capture System. These instructions for use must be kept in a convenient location providing direct and ready access for the Operating Room and biomedical technical personnel at any time.**

WARNING: **CONTRAfuran / SENSOfluran CONNECT are not designed to adsorb, capture or detect nitrous oxide. Since **CONTRAfuran** / SENSOfluran CONNECT CANNOT adsorb nitrous oxide, the system must be connected correctly to the Operating Room AGSS (Anesthetic Gas Scavenging System) to ensure safe and complete removal of nitrous oxide anesthetic gas.**

WARNING! Follow strictly these instructions for use

Table of Contents

1. What is **CONTRAfuran** Anesthetic gas canister and what is the SENSOfluran CONNECT Fill level control unit?
2. Appropriate use
3. Safety instructions for the **CONTRAfuran** / SENSOfluran CONNECT System
4. Transportation, storage, and operation of **CONTRAfuran** Anesthetic gas canister and the SENSOfluran CONNECT Fill level control unit
 - 4.1. Transportation
 - 4.2. Storage of new (unused) **CONTRAfuran** Anesthetic gas canisters
 - 4.3. Operation
5. Application area of the **CONTRAfuran** / SENSOfluran CONNECT Anesthetic Gas Capture System in surgical Operating Rooms with Anesthetic Gas Scavenging System (AGSS) exhaust facilities
6. When and how to replace the used **CONTRAfuran** Anesthetic gas canister
7. Return and handling of the used (filled) **CONTRAfuran** Anesthetic gas canister
8. Maintenance and calibration interval of the SENSOfluran CONNECT Fill level control unit
9. Technical data
10. Customer service

1. What is CONTRAfluran Anesthetic gas canister and what is the SENSOfluran CONNECT Fill level control unit?

CONTRAfluran / SENSOfluran CONNECT System is a new, innovative technology protected by international patents and patent applications with which fluorinated halocarbon inhaled anesthetic gases such as sevoflurane, desflurane and isoflurane administered during surgery are completely adsorbed and removed from the exhaust of the anesthetic machine. The CONTRAfluran Anesthetic gas filter canister contains a unique adsorber material characterized by its highly specific and controlled physical and chemical properties such as grain structure, porosity and surface area. The adsorber is highly selective and efficient in retaining volatile fluorinated halocarbon anesthetic gases which enter the canister. The SENSOfluran CONNECT Fill level control unit functions as a canister holder and is equipped with a sensor to monitor the fill level (adsorbed gas capacity) of the CONTRAfluran Anesthetic gas canister. SENSOfluran CONNECT indicates when the CONTRAfluran Anesthetic gas canister is filled and must be changed.

WARNING: CONTRAfluran / SENSOfluran CONNECT System CANNOT adsorb, capture or detect nitrous oxide. Since CONTRAfluran CANNOT adsorb nitrous oxide, the system must be connected correctly to the Operating Room AGSS gas removal system to ensure safe removal of non-adsorbed gases. CONTRAfluran and SENSOfluran CONNECT are only indicated for with the adsorption of halogenated anesthetic gases (sevoflurane, desflurane and isoflurane).

2. Appropriate use

The CONTRAfluran Anesthetic gas canister is intended to adsorb and remove halogenated anesthesia gases (such as sevoflurane, desflurane and isoflurane). The CONTRAfluran / SENSOfluran CONNECT Anesthetic Gas Capture System is designed exclusively for use with AGSS-equipped Operating SENSOfluran environments. The SENSOfluran CONNECT Fill level control unit must only be used in combination with the CONTRAfluran Anesthetic gas canister during induction and maintenance of inhaled anesthesia in patients using volatile anesthetics.

CONTRAfluran and SENSOfluran CONNECT may safely be used in situations where N₂O (nitrous oxide) sedation is being administered however, they are not intended to detect or capture nitrous oxide or other non-halogenated anesthetic gases!

Please note that the CONTRAfluran Anesthetic gas canister and the SENSOfluran CONNECT Fill level control unit must only be installed and used as instructed by qualified personnel.

Strict compliance with the directions included in these Instructions for use will ensure the safe and efficient use of the CONTRAfluran / SENSOfluran CONNECT System. Only the ZeoSys-supplied hoses or other identical medical connectors should be used, forming a tight seal to the 22mm tapered connector inlet port. The tight seal shall be installed and permanently fixed by suitable tools. ZeoSys GmbH or its authorized representatives will not accept any liability or damage, which is caused by incorrect or unreasonable use of the CONTRAfluran Anesthetic gas canister or the SENSOfluran CONNECT Fill level control unit.

3. Safety instructions for the CONTRAfluran / SENSOfluran CONNECT System

The following safety instructions must always be followed during the operation, maintenance and repair work performed on this device.

Failure to comply with these safety instructions can lead to danger for the operating personnel, and damage to the device!

- The SENSOfluran CONNECT Fill level control unit must never be opened and no modifications should be made to it. The exchange of components and repairs must be done by qualified / trained personnel.
- The CONTRAfluran / SENSOfluran CONNECT System is to be operated within its specified temperature range from +5°C to +35°C.
- The SENSOfluran CONNECT Fill level control unit must be connected to and powered using only the provided plug-in power supply unit. It is designed for an operating voltage of 100V -240V AC / 47 -63 Hz.

- Caution should be exercised when installing and operating the **CONTRAfluran / SENSOfluran CONNECT System** in an area where there is significant exposure to explosive and / or combustible materials.
- Damaged or inoperable or defective part of the **CONTRAfluran / SENSOfluran CONNECT System** or the entire medical device must be removed from operation and secured from unauthorized use. Return all such systems or devices to ZeoSys Medical immediately for replacement.
- **SENSOfluran CONNECT Fill level control unit** and **CONTRAfluran Anesthetic gas canister** contain no user device serviceable parts or components and should be returned immediately to ZeoSys Medical for repair and/ or replacement.
- It must be ensured that the medical inlet and outlet hose connections of the **CONTRAfluran / SENSOfluran CONNECT System** provide a tight seal, and that the anesthesia workstation/ machine is configured to operate in **passive scavenging mode**.
- The examples of anesthesia machines listed below do not currently allow this passive mode, or allow it only through device intervention, and are therefore **NOT** suitable for the intended use of the **CONTRAfluran Anesthetic gas canister** with the **SENSOfluran CONNECT Fill level control unit**.

-Dräger Atlan

-Getinge anesthesia machines

-Dameca Philips anesthesia machines

- For the following anesthesia machines, a conversion into passive mode is only possible with the installation of an additional part:

-Dräger Perseus

-GE (Avance, Aisys, Aestiva, Aespire, Carestation 600/650)

-Mindray A8/A9

In this case the representative of the respective anesthesia machine company must be consulted for installation of the additional part needed for converting the specific anesthesia machine model into passive mode.

In case of questions, please contact the anesthesia machine manufacturer.

- Upon connection of **CONTRAfluran / SENSOfluran CONNECT System**, the anesthesia workstation should be put through the manufacturers recommended start-up / operational testing (e.g. "self-test"). Any test errors should be corrected and evaluated to confirm that they do not result from the installation of **CONTRAfluran / SENSOfluran CONNECT System**, prior to proceeding with the use of the workstation.

4. Transportation, storage and operation of **CONTRAfluran Anesthetic gas canister and **SENSOfluran CONNECT Fill level control unit****

This section describes the transportation, storage, and operation of the **CONTRAfluran Anesthetic gas canister** and the **SENSOfluran CONNECT Fill level control unit**.

To ensure the performance of your **CONTRAfluran Anesthetic gas canister** and **SENSOfluran CONNECT Fill level control unit** and in order to help ensure your safety, be sure to follow all recommended application procedures.

4.1. Transportation

During the transportation of the **CONTRAluran / SENSOfluran CONNECT System**, the following points should be considered:

- The shipping carton should be transported upright with the label showing right-side up.
- The shipping carton should not be inverted or stored on its side.
- Carefully open the carton from the top avoiding damage to the contents.
- Retain the original **CONTRAluran** shipping carton, red caps and zipper bags for sealing and return transportation of filled canisters.

4.2 Storage of new (unused) **CONTRAluran Anesthetic gas canisters**

For storage of new **CONTRAluran** Anesthetic gas canisters and / or **SENSOfluran CONNECT** Fill level control units which are not currently in use, please consider the following:

- The shipping carton should be stored with the label showing right-side up.
- Packaging materials including shipping carton, zipper bags and red caps should be retained for the return of canisters.
- Store in a dry, clean location at room temperature.

4.3 Operation

Before using the **CONTRAluran / SENSOfluran CONNECT System**, please ensure that all components needed for proper installation and use are available.

This includes:

- The **CONTRAluran** Anesthetic gas canister
- The **SENSOfluran CONNECT** Fill level control unit
- Accessories: flexible anesthesia hose (22mm ISO taper connector) and operational adapters for use with different exhaust ports
- Mounting rail splint at the anesthesia machine for attaching the **CONTRAluran / SENSOfluran CONNECT System**

Prior to use, confirm that the anesthesia machine has been converted to **passive scavenging mode**.

A1. *Securely fasten the **SENSOfluran CONNECT** Fill level control unit to the mounting rail of the anesthesia workstation / machine.*

A2. Connect the mains adapter into a suitable power outlet.

Note: Be sure that the mains adapter is connected to an easily accessible power outlet. Unplug the mains adapter from the power outlet in case of emergency.

The SENSOfluran CONNECT Fill level control unit enters a self-diagnostic mode and performs a self-test consisting of the following steps:

- 1 second visual and acoustic test
 - the green and yellow LEDs will illuminate
 - an acoustic signal sounds
- 10 second test of the optical presence sensor for the canister
- 5 minutes warm-up period
 - the green LED flashes

At the end of this self-test, 3 LEDs light up one after the other and an acoustic signal sounds.

(This means that no **CONTRAfuran** Anesthetic gas canister has been placed into the SENSOfluran CONNECT Fill level control unit)

Follow the next steps to continue and complete the connection of the **CONTRAfuran** / SENSOfluran CONNECT System to the anesthesia machine as shown schematically on the following pictures:

A3. Place a **CONTRAfuran** Anesthetic gas canister with an intact red safety seal into the SENSOfluran CONNECT Fill level control unit. Break the safety seal by removing (and retaining) the red seal cap from the **CONTRAfuran** Anesthetic gas canister. The green (bottom) LED flashes, indicating that the sensor is in warm-up period of 15 minutes. Wait till this warm-up period is complete. The bottom LED will return to a steady green light.

A4. Connect the **CONTRAfuran** Anesthetic gas canister inlet to the AGS outlet (exhaust gas outlet) connector of the anesthesia workstation using the flexible hose provided (22 mm ISO taper connector).

A5. Connect the Visual Flow Indicator outlet port of the SENSOfluran CONNECT Fill level control unit to the AGSS wall or ceiling connection in the Operating Room as shown in Section 5. Make sure that the connection is tight and secure. Check the operational readiness of AGSS (Anesthetic Gas Scavenging) connection with the help of the Visual Flow Indicator.

A6. Record the anesthetic gas(es) being collected by checking the appropriate box on the top of the canister.

B. Possible LED Display Signals of SENSOfluran CONNECT Fill level control unit

B1. Fill level indication

The differently colored LEDs (green, yellow, yellow) of the device indicate the quality of the filtered outlet gas and thus the fill level of the canister.

NOTE: The top-most LED of 4 LEDs (red) does not illuminate under any conditions of normal use.

GREEN LED (steady ON, bottom)

Indicates that the CONTRAfluran Anesthetic gas canister is adsorbing halogenated anesthesia gases exhausted from anesthesia work station (other gases pass through and are removed by the AGSS).

GREEN/YELLOW LED (bottom)

The filter is still functioning and filtering out anesthetic gas, however the CONTRAfluran Anesthetic gas canister is nearing its capacity.

YELLOW LED (bottom)

The capacity of the filter is almost exhausted. The canister change must be prepared. Provide a new the CONTRAfluran Anesthetic gas canister for canister change.


YELLOW flashing LED + audible interval signal (top)

The capacity of the filter is exhausted.

The CONTRAfluran Anesthetic gas canister must immediately be replaced with a new one.


B2. Error messages:

GREEN / YELLOW / YELLOW

3 LEDs light up one after the other and an acoustic signal sounds.

This means that no canister has been placed in the SENSOfluran CONNECT Fill level control unit. Place the **CONTRAfuran** Anesthetic gas canister in the SENSOfluran CONNECT Fill level control unit and wait for the green continuous light.

FLASHING GREEN

After a canister change, the green (bottom) LED flashes, indicating that the sensor is in the warm-up period. When the warm-up period of 15 minutes is complete, the bottom LED will return to a steady green light.


NOTE: The **CONTRAfuran** Anesthetic gas canister remains operational during the warm-up period.

GREEN FLASHING WITH ACOUSTIC SIGNAL

If the green LED flashes and an acoustic signal sounds at the same time, the gas sensor is defective. Please remove from service and arrange replacement with the customer service.

5. Application area of the **CONTRAfuran / SENSOfluran CONNECT Anesthetic Gas Capture System in surgical Operating Rooms with Anesthetic Gas Scavenging System (AGSS) exhaust facilities.**

6. When and how to replace the used CONTRAfluran Anesthetic gas canister

YELLOW LED (bottom)

The capacity of the filter is almost exhausted. The canister change must be prepared. Provide a new CONTRAfluran Anesthetic gas canister for canister change.

YELLOW flashing LED + audible interval signal (top)

The capacity of the filter is exhausted.

The CONTRAfluran Anesthetic gas canister must immediately be replaced with a new one.

Place a CONTRAfluran Anesthetic gas canister with an intact red safety seal into the SENSOfluran CONNECT Fill level control unit. Break the safety seal by removing (and retaining) the red seal cap from the CONTRAfluran Anesthetic gas canister.

The green (bottom) LED flashes, indicating that the sensor is in a warm-up period of 15 minutes.

After completion of this warm-up period the bottom LED will return to steady green light.

7. Return and handling of the used (filled) CONTRAfluran Anesthetic gas canister

Where possible, replace the red cap onto the used CONTRAfluran Anesthetic gas canister.

Prior to storage and return, please ensure that the exterior surface of the CONTRAfluran Anesthetic gas canister is wiped with a cloth dampened with an aqueous-based cleaning solution (ZeoSys Medical recommends the use of hydrogen peroxide wipes, concentration 1% -1.5%).

DO NOT USE alcohol or aldehyde-based cleaning agents as they will damage the SENSOfluran CONNECT sensor.

- The canister must then be placed into the protective zipper bag and **securely sealed**.
- Six used canisters, securely sealed in their zipper bags, should be packed into the **original shipping carton** and stored in a cool, dry, well-ventilated place, pending return to ZeoSys Medical.
- Before returning the filled used canisters to ZeoSys Medical please label the **shipping carton** with the name and address of your hospital.

8. Maintenance and calibration interval of SENSOfluran CONNECT Fill level control unit

Follow hospital-established protocol for cleaning the anesthesia workstation to determine cleaning frequency for SENSOfluran CONNECT Fill level control unit.

To clean the surface of the SENSOfluran CONNECT Fill level control unit, wipe down with a cloth dampened with hydrogen peroxide (1%-1.5%) or other aqueous biocidal detergent. **DO NOT USE** alcohol or aldehyde-based cleaners, as they will damage the sensor and negatively impact the performance of the SENSOfluran CONNECT Fill level control unit.

CAUTION: Avoid excessive „wetting“ near the electrical components or the sensor during cleaning period or sanitization activities! Use non-alcohol-based wipes only!

Allow at least 10 minutes for exposure and drying time before placing the SENSOfluran CONNECT Fill level control unit back into operation.

The gas sensor unit within the SENSOfluran CONNECT Fill level control unit must be re-calibrated every **12 months**.

Calibration can only be conducted by ZeoSys Medical GmbH. Please refer to the **calibration sticker** showing the re-calibration date applied to the side of SENSOfluran CONNECT Fill level control unit. For a newly-calibrated unit, please contact Customer Service at least 1 month prior to expiration date.

9. Technical data

CONTRAfuran	Product-class according to EU MDR	Class 1
Anesthetic gas canister	Temperature	Operation: +5°C to +35°C Storage: -5°C to +35
	Relative humidity	Operation: to 70% Storage: to 70%
	Storage capacity	Approx. 400g
	Flow-Resistance	≤ 0,15 mbar
	Height	19 cm
	Diameter	12 cm
	Weight	Approx. 1000 g
	Volume	2 liters
SENSOfuran	External material	Powder-coated aluminum
CONNECT	Temperature	+5°C to +35°C
Fill level control unit	Voltage	Power supply 100V-240V AC / 47-63Hz Output 6,0V DC Power consumption approx: 6W
	Weight	Approx. 1400g
	Gas sensor circuit board	(50x55x20) mm
	Circuit board coating	Bectron MR 3404
	External wall mains power adapter	GSM06E06 -P1J; GEM06 I6-P1J

Products	Part numbers
SENSOfuran CONNECT EU	Zeo 000073
SENSOfuran Starter Kit EU	Zeo 000173
SENSOfuran CONNECT UK	Zeo 000074
SENSOfuran Starter Kit UK	Zeo 000174
SENSOfuran CONNECT US	Zeo 000075
SENSOfuran Starter Kit US	Zeo 000175
Electronic board	Zeo 000042
CONTRAfuran™ Anesthetic gas canister	Zeo 000050
Flexible hose 22mm (interface 40 cm, length 3,20m)	1574000
Connector 22M-30F	1971000
Connector 22M-30M	1970000
Connector 22M-22F	1961000
Flow indicator	Zeo 000072
Power supply unit	Zeo 000053
Power adapter EU	Zeo 000054
Power adapter UK	Zeo 000055
Power adapter US	Zeo 000056

Manufacturer's declaration

Applicable Classifications per IEC 60601-1, Clause 6:

Classification	
Protection against electric shock	Class II
Suitability for an oxygen rich environment (greater than 25 % for ambient pressures up to 110 kPa)	Suitable
Mode of operation	Continuous

EMC Compliance

The **SENOfluran CONNECT Fill level control unit** complies with IEC 60601-1-2:2014+AMD1:2020. However, special precautions concerning EMC must be taken for all medical electrical equipment.

WARNING: Use of the **SENOfluran CONNECT Fill level control unit** adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the **SENOfluran CONNECT Fill level control unit** and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **SENOfluran CONNECT Fill level control unit**, including cables specified by Baxter. Otherwise, degradation of the performance of this equipment could result.

NOTE: The **SENOfluran CONNECT Fill level control unit** may display inconsistent indications of the filter fill level in the presence of EM disturbances (e.g. high-frequency surgical equipment). Once these EM disturbances stop the **SENOfluran CONNECT Fill level control unit** will perform as intended.

Emissions and Immunity

Table 1 - Emission limits and compliance

Phenomenon	Emission Limits	Compliance
Conducted and radiated RF EMISSIONS	CISPR 11	Group 1, Class B

Table 2 - Immunity test levels (Enclosure Port)

Phenomenon	Basic EMC standard or test method	Immunity Test Levels
Electrostatic Discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 60 Hz
Proximity magnetic fields	IEC 61000-4-39	See table 4

Table 3 - Test levels -proximity fields from RF wireless communication equipment

Test Frequency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Immunity Level (V/m)
385	380–390	TETRA 400	Pulse 18 Hz	1.8	0.3	27
450	430–470	GMRS 460 FRS 460	Frequency ±5 kHz 1 kHz Sine	2	0.3	28
710 745 780	704–787	LTE Band 13,17	Pulse 217 Hz	0.2	0.3	9
810 870 930	800–960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse 18 Hz	2	0.3	28
1720 1845 1970	1700–1990	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1,3, 4, 25 UMTS	Pulse 217 Hz	2	0.3	28
2450	2400–2570	Bluetooth WLAN 802.11 (b/g/n) RFID 2450 LTE Band 7	Pulse 217 Hz	2	0.3	28
5240 5500 5785	5100–5800	WLAN 802.11 (a/n)	Pulse 217 Hz	0.2	0.3	9

Table 4 –Text levels - Proximity magnetic fields

Test Frequency	Modulation	Immunity Level (A/m)
30 kHz	CW	8
134.2 kHz	Pulse 2.1 kHz	65
13.56 MHz	Pulse 50 kHz	7.5

Table 5 - Immunity test levels (a.c. port)

Phenomenon	Basic EMC standard or test method	Immunity Test Levels
Electrical fast transients / bursts	IEC 61000-4-4	$\pm 2 \text{ kV}$ 100 kHz repetition frequency
Surges Line-to-line	IEC 61000-4-5	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}$
Surges Line-to-ground	IEC 61000-4-5	$\pm 0.5 \text{ kV}, \pm 1 \text{ kV}, \pm 2 \text{ kV}$
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips	IEC 61000-4-11	0 % U_t ; 0.5 cycle At $0^\circ, 45^\circ, 90^\circ, 135^\circ, 180^\circ, 225^\circ, 270^\circ$ and 315° 0 % U_t ; 1 cycle and 70 % U_t ; 25/30 cycles Single phase; at 0°
Voltage interruptions	IEC 61000-4-11	0 % U_t ; 250/300 cycle

Note: UT is the a.c. mains voltage prior to application of the test level.

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